

INTRODUCTION. Susceptibility testing, identification by DNA gene sequencing and DNA fingerprinting of the rapidly growing mycobacteria and other nontuberculous mycobacteria (NTM) and related aerobic actinomycetes are performed at the UTHSCT Mycobacteria/Nocardia Laboratory. The laboratory is College of American Pathologists (CAP) accredited and holds a CLIA license, a Florida CLIA license, and a Pennsylvania license.

The Mycobacteria/Nocardia Laboratory has been in operation for approximately 40 years and accepts pure culture isolates (not gross specimens) of the above groups for susceptibility, identification, and DNA fingerprinting. **Please note there will be an additional charge and increased turn around time for any isolates that are not pure.** If desired, gross specimens should be submitted to the UTHSCT Pathology Laboratory (See Mycobacterium/Mycology Referral; call **903-877-5745** for details).

SUSCEPTIBILITY TESTING. Susceptibility testing is performed using the CLSI recommended broth microdilution MIC method for the NTM, Nocardia, and related aerobic actinomycetes. For rapidly growing mycobacteria (RGM), the routine panel of drugs* includes clarithromycin, amikacin, imipenem, linezolid, tobramycin (for *M. chelonae/immunogenum* complex), cefoxitin, moxifloxacin, ciprofloxacin, trimethoprim-sulfamethoxazole, and doxycycline and/or minocycline. Incubation time is 3-5 days at 30°C except for clarithromycin which is held up to 14 days to check for inducible in vitro resistance (*erm* gene). However, if isolates of *Mycobacterium abscessus* complex, *M. chelonae*, *M. immunogenum* and *M. mucogenicum* complex are sequenced in our laboratory, extended incubation will not be necessary and the clarithromycin MICs will be reported after 3-4 days incubation as susceptible or resistant based on sequencing of the *erm* gene, *rpoB* sequence identification, and/or exclusion of a 23S rRNA gene mutation.

In isolates with resistant clarithromycin MICs, macrolide mutational resistance can be confirmed by ordering macrolide mutational sequencing 23S rRNA (Test Code 13). Additionally, confirmation of amikacin mutational resistance can be performed in isolates with high level amikacin MICs by ordering 16S rRNA mutational sequencing (Test Code 13).

For Nocardia and other related aerobic actinomycetes, the panel of drugs* includes amoxicillin/clavulanic acid, linezolid, clarithromycin, trimethoprim-sulfamethoxazole, ciprofloxacin, moxifloxacin, ceftriaxone, imipenem, tobramycin, amikacin, doxycycline and/or minocycline. Incubation time is 3-5 days at 35°C for most species.

Slowly growing NTM (other than *M. avium* complex and *M. kansasii*) are tested against clarithromycin, ethambutol, amikacin, trimethoprim-sulfamethoxazole, rifabutin, ciprofloxacin, moxifloxacin, minocycline, linezolid, and rifampin.

Current recommendations by the ATS and the CLSI/NCCLS M24-A2, 2011 advise testing only clarithromycin and amikacin (proposal made to report this latter antibiotic) against isolates of *Mycobacterium avium* complex (MAC) as these are the only agents which have been shown to have *in vitro* correlation of MICs with clinical response. Clarithromycin is used as a "class drug" in susceptibility testing of azithromycin and other related macrolides. Thus, isolates of MAC are tested only for susceptibility to clarithromycin and amikacin. Resistance to clarithromycin confers resistance to azithromycin and vice versa. It may be reasonable to test agents such as moxifloxacin and linezolid and these may, in some cases, provide useful adjunctive treatment options. They should not, however, be used as treatment substitutions for any of the standard treatment agents (macrolide, ethambutol, rifampin, rifabutin) and their efficacy for treatment remains unproven. First line TB agents (ethambutol, rifampin, and isoniazid) are not reported with isolates of MAC in accordance with CLSI and ATS recommendations. Currently, there are no CLSI recommended breakpoints for amikacin with MAC, however, a proposal has been made to consider amikacin MICs ≥ 64 $\mu\text{g/mL}$ (IV treatment) and > 64 $\mu\text{g/mL}$ (inhaled treatment) as resistant so that resistance would be defined as ≥ 64 $\mu\text{g/mL}$ depending upon method of administration of amikacin. Mutational resistance can be confirmed for clarithromycin and/or amikacin by ordering Test Code 13 and specifying amikacin, clarithromycin or both.

In accordance with the American Thoracic Society (ATS) and CLSI recommendations, isolates of *M. kansasii* are reported with rifampin and clarithromycin susceptibility only, if they are rifampin susceptible. Additional agents can be included with special physician requests. Rifampin resistant *M. kansasii* will include a panel of drugs (same as noted for slowly growing NTM other than MAC). Incubation time is 7-14 days at 35°C.

The methods and breakpoints for susceptibility testing of RGM and aerobic actinomycetes have been published by the CLSI. Recommendations for breakpoints for these organisms using the broth microdilution MIC method with nine drugs (amikacin, tobramycin, trimethoprim-sulfamethoxazole, cefoxitin, imipenem, linezolid, doxycycline/minocycline, clarithromycin and ciprofloxacin) are shown in Tables 1,2, and 4. Breakpoints for tigecycline have not been addressed by the CLSI yet for NTM or other aerobic actinomycetes. However, our laboratory will report MICs to tigecycline for RGM without interpretation upon special request for an additional charge.

Generally, turn-around time for susceptibility testing is approximately **14-30**

working days for RGM (includes macrolide induction testing), **7 – 14 working days** for *Nocardia* and related aerobic actinomycetes, and **21-35 working days** for the slowly growing NTM and some other aerobic actinomycetes.

Please note, if isolates are identified in our laboratory as containing a functional *erm* gene (i.e. most *M. abscessus* subsp. *abscessus*, *M. abscessus* subsp. *bolletii*, these isolates will be reported as macrolide resistant and 14d clarithromycin induction will not be required. T-A-T for susceptibility will be reduced to approximately 7-10 working days, provided the culture submitted is a pure culture isolate and no problems occur with growth (Nash, et al., AAC 2006;50:3476; Nash, et al., JAC 2005;55:170; Brown-Elliott et al., 2015;53:1211).

If the isolate is identified in our laboratory as *M. chelonae*, *M. abscessus* subsp. *massiliense*, *M. immunogenum*, *M. mucogenicum* complex or a sequence variant of *M. abscessus* which does not harbor a functional *erm* gene, it will be reported as macrolide susceptible without 14d induction (Brown-Elliott et al., JCM 2015;53:875). Again, T-A-T will be reduced to approximately 7-10 working days with the same provisions as described above.

Isolates of *Mycobacterium haemophilum* require longer incubation with an average turn-around-time of approximately 6 weeks. The susceptibility panel includes amikacin, ciprofloxacin, clarithromycin, doxycycline and/or minocycline, linezolid, rifampin, and TMP-SMX. The testing uses the CLSI recommended agar disk elution method with hemin supplementation (CLSI, 2011).

The turn-around-times given assume receipt of a viable pure culture isolate. Results are sent by FAX and mail. **FOR SUSCEPTIBILITY RESULTS OR TO CHECK RECEIPT OF ISOLATES PLEASE CALL (903) 877-7978.**

Recent studies have shown that previously treated isolates of *M. abscessus* may be difficult to grow adequately for susceptibility testing in broth, probably due to antibiotic pressure and stress on the isolates. Extended incubation, which may be deleterious to some antimicrobials may be required and if so, this will be noted on the report. If the organism does not grow within 5 days, the report will state that susceptibility testing is unable to be performed within the acceptable time period and that molecular sequencing of the 16S rRNA gene and *erm* gene are recommended to determine amikacin and clarithromycin susceptibility/resistance respectively. (Amikacin and clarithromycin are two important agents used for treatment of *M. abscessus* complex.) Clinical consultation for treatment recommendations and interpretation of these results should be considered.

If the 14 day clarithromycin MIC result is 4 µg/mL (intermediate), the submitter should consider ordering *erm* gene sequencing to more accurately determine susceptibility/resistance of the isolate.

More studies are needed to assess the utility of sequencing of the *erm* gene among isolates of the *M. fortuitum* group except for isolates of *M. senegalense* and *M. peregrinum* which do not contain functional *erm* genes. Until such evaluations are performed, extended (up to 14 day) incubation of isolates to assess macrolide susceptibility/ resistance will be required unless the isolate has an initial (3-5 day) MIC of ≥ 8 $\mu\text{g}/\text{mL}$ (suggesting mutational resistance which can be confirmed by ordering Test code 13 for clarithromycin).

IDENTIFICATION. The application of molecular techniques for the identification of mycobacteria and nocardia has become the “gold standard” among methods of identification for these species.

The Mycobacteria/Nocardia Laboratory primarily uses *rpoB* sequencing along with *erm* gene sequence for identification of most clinically significant non-pigmented rapidly growing mycobacteria (Adékambi 2003; Steingrube 1995, 1997; Nash 2009). However, for slowly growing mycobacteria, Nocardia and other aerobic actinomycetes, sequence analysis of 500 bp segment of the 16S rRNA gene is also often used. Sequencing of the *sec A* gene is also used for the identification of Nocardia (Conville JCM 2006; Brown-Elliott CMR 2006). Sequencing of the entire 16S rRNA gene (approximately 1500 bp) is not a practical method for routine test identification and thus sequencing of other target genes may be necessary. If species identification is not possible using the 16S rRNA, our laboratory employs other target genes (e.g. *rpoB*, *sec A*, etc.) or in some cases full 16S rRNA gene sequencing. The interpretation of gene sequences follows the recommendations of the CLSI published in the MM18-A (CLSI, 2008). Isolates which give indeterminate results may be submitted for additional testing which may include full 16S rRNA sequence and/or multi-gene sequence analysis if required.

Recent studies have shown interspecies gene transfer among species and subspecies of RGM especially within the *M. abscessus/chelonae* complex and *M. fortuitum* group. This new finding necessitates the use of multiple gene sequences in order to determine definitive identification of species level. This additional step adds a significant amount of work to the identification process and we will be re-assessing our turn-around-times (based on working days) as time to final test results will be longer. Preliminary test results should, however, be available in approximately 72 hours. Please call us if you have an urgent need for isolate ID and we will attempt to identify as soon as possible during this transition period.

FOR SEQUENCING OR DNA STRAIN TYPING RESULTS, PLEASE CALL (903) 877-5947 or (903) 877-7683.

DNA FINGERPRINTING. Isolates of NTM or other aerobic actinomycetes may also be submitted for DNA fingerprinting by pulsed field gel electrophoresis (PFGE) in cases of outbreaks, pseudo-outbreaks, or epidemics. Variable number tandem repeat (VNTR) is also available for strain typing of *Mycobacterium avium* and *M. intracellulare* (Iakhiaeva et al., 2013; Iakhiaeva et al., 2016). Before sending isolates for restriction fragment length polymorphism (RFLP) analysis, please call for consultation. Turn-around-time for PFGE is approximately 6-8 weeks and 1-2 weeks for VNTR from receipt of the isolate depending on the type and number of organisms submitted for testing. (See [Fees for Laboratory Services](#))

FOR TECHNICAL CONSULTATION CALL BARBARA BROWN-ELLIOTT AT (903) 877-7685.

***Subject to Change (See Table 3 for expanded current panels available)**

Table 1. Suggested broth microdilution breakpoints for rapidly growing mycobacteria^a.

Antimicrobial Agent	Minimal Inhibitory Concentration (µg/mL) for Category		
	Susceptible	Intermediate	Resistant
Amikacin ^b	≤16	32	≥64
Cefoxitin	≤16	32-64	≥128
Ciprofloxacin	≤1	2	≥4
Clarithromycin	≤2	4	≥8
Doxycycline/Minocycline	≤1	2-4	≥8
Imipenem/Meropenem	≤4	8-16	≥32
Linezolid	≤8	16	≥32
Moxifloxacin ^c	≤1	2	≥4
Tobramycin ^d	2	4	≥8
Trimethoprim- Sulfamethoxazole	≤2/38	-	≥4/76

^aBreakpoints from the CLSI (formerly NCCLS), NCCLS-M24-A, 2003; CLSI M24-A2, 2011.

^bAmikacin resistance ≥64 µg/mL applies to IV treatment; >64 µg/mL is considered resistant for inhaled treatment. (Amikacin resistance can also be confirmed by sequencing the 16S rRNA gene for the mutation, Test Code 13).

^cMoxifloxacin breakpoints tentative in M24-2, 2011.

^dTobramycin is only reported for *M. chelonae*/*M. immunogenum*.

NOTE: There are no interpretive criteria for tigecycline established with nontuberculous mycobacteria currently. This antibiotic is tested upon special request. MICs are given without interpretation of values.

Table 2. CLSI/NCCLS M24-A2, Suggestions for susceptibility testing of *M. abscessus*, *M. chelonae*, and the *M. fortuitum* group by broth microdilution^a.

Drug	Comment
Tobramycin	If the initial MIC is >4 µg/mL, the test should be repeated. If the repeat result is >4 µg/mL, the MIC should be reported with a comment. ^b
Sulfonamides	MIC is read at 80% inhibition of growth.
Doxycycline	Breakpoints are ≤1 µg/mL (susceptible), 2-4 µg/mL (intermediate), and ≥8 µg/mL (resistant).
Cefoxitin	Breakpoints are ≤16 µg/mL (susceptible), 32-64 µg/mL (intermediate), and ≥128 µg/mL (resistant).
Imipenem	If MIC for <i>M. fortuitum</i> group, <i>M. smegmatis</i> group, or <i>M. mucogenicum</i> group is >8 µg/mL, test should be repeated with incubation period of no more than 3 days. If the repeat result is >8 µg/mL, the MIC should be reported with comment. ^b
Amikacin	Isolates of <i>M. abscessus</i> for which MIC is ≥64 µg/mL should be retested. If the repeat result is ≥64 µg/mL, the MIC should be reported with a comment. ^b
Clarithromycin	Isolates of <i>M. fortuitum</i> group with a trailing endpoint should be considered resistant. Extended incubation up to 14 days should be performed to detect the inducible <i>erm</i> gene.

^a For laboratories that infrequently isolate rapidly growing mycobacteria, sending isolates to an experienced reference laboratory is recommended. For laboratories that perform MIC testing, (i) proficiency testing by comparison of test results with those of an experienced reference laboratory is necessary upon initial validation and at regular intervals thereafter and (ii) identification of isolates to the species level or, at a minimum, differentiation of the *M. fortuitum* group from the *M. chelonae* and *M. abscessus* group is recommended.

^b Comment: (i) the MIC is greater than expected for this species and (ii) if the drug is being considered for therapy, the laboratory should be notified so the isolate can be sent to a reference laboratory for confirmation of resistance. Please note that isolates of *M. immunogenum* are usually resistant to tobramycin (MIC >4 µg/mL) in contrast to isolates of *M. chelonae* which usually have tobramycin MICs ≤2 µg/mL. Therefore, sequence identification may also be helpful.

Table 3. Antibiotic Panels Available

Slowly Growing Mycobacteria	Rapidly Growing Mycobacteria	Nocardia
Amikacin	Amikacin	Amikacin
Ciprofloxacin	Cefoxitin	Amoxicillin-Clavulanic Acid (Augmentin)
Clarithromycin ¹	Ciprofloxacin	Ceftriaxone
Doxycycline or Minocycline	Clarithromycin	Ciprofloxacin
Ethambutol ²	Doxycycline	Doxycycline
Linezolid ^{3,6}	Clarithromycin	Imipenem
Moxifloxacin ⁴		Linezolid
Rifabutin	Imipenem	Minocycline
Rifampin	Linezolid	Moxifloxacin ⁴
TMP/SMX	Minocycline	TMP-SMX
	Moxifloxacin ⁴	Tobramycin
	TMP-SMX	
	Tigecycline ⁵	
	Tobramycin ⁶	

¹ Class drug for newer macrolides (i.e., azithromycin and clarithromycin)

² Please note that the MIC testing of ethambutol (which is considered to be technically difficult and unreliable by many investigators) is currently being discussed by the CLSI and this agent may be removed dependent upon the pending CLSI consensus report.

³ MICs for linezolid interpreted using RGM breakpoints until final breakpoints have been addressed for slowly growing species.

⁴ MIC interpretations based on CLSI M100-S15, 2005, bacterial breakpoints and tentative in M24A2.

⁵ MIC reported without interpretation since no breakpoints yet established. Testing performed upon special request.

⁶ MIC only reported for *M. chelonae/immunogenum* complex.

*Isolates of *Mycobacterium avium* complex (MAC) are routinely tested for clarithromycin and amikacin susceptibility only, and rifampin susceptible *M. kansasii* are tested for susceptibility to rifampin and clarithromycin only.

**Upon physician request, isolates of rifampin susceptible *M. kansasii* may be tested for susceptibility to other agents.

⁶ Upon request, isolates of MAC may be tested for susceptibility to moxifloxacin and linezolid but no first line TB drugs are reported.

NOTE: Drugs tested are subject to change without notice!

For availability of susceptibilities to agents other than listed, please call (903) 877-7685.

Table 4. Suggested broth microdilution interpretive criteria for *Nocardia* and other aerobic actinomycetes¹.

Antimicrobial Agent	Minimal Inhibitory Concentration ($\mu\text{g/mL}$) for Category		
	Susceptible	Intermediate	Resistant
Amikacin	≤ 8	-	≥ 16
Amoxicillin/clavulanic acid	$\leq 8/4$	16/8	$\geq 32/16$
Ceftriaxone	≤ 8	16-32	≥ 64
Ciprofloxacin	≤ 1	2	≥ 4
Clarithromycin	≤ 2	4	≥ 8
Imipenem	≤ 4	8	≥ 16
Linezolid	≤ 8	-	-
Minocycline/Doxycycline	≤ 1	2-4	≥ 8
Moxifloxacin	≤ 1	2	≥ 4
Trimethoprim- Sulfamethoxazole	$\leq 2/38$	-	$\geq 4/76$
Tobramycin	≤ 4	8	≥ 16

¹Reference for tentative breakpoints: CLSI, M24-A-2, 2011.
For isolates of *R. equi*, use the CLSI, M100 criteria for *S. aureus* and include results for vancomycin and rifampin.

Table 5. Suggested broth microdilution breakpoints for slowly growing nontuberculous mycobacteria.¹

Antimicrobial Agent	Minimal Inhibitory Concentration ($\mu\text{g/mL}$) for Category		
	Susceptible	Intermediate	Resistant
Amikacin ^{2,3}	≤ 16	32	≥ 64
Rifampin	≤ 1	-	≥ 2
Ethambutol	≤ 2	4	≥ 8
Rifabutin	≤ 2	-	≥ 4
Ciprofloxacin	≤ 1	2	≥ 4
Sulfamethoxazole/TMP-SMX	≤ 32 (2/38)	-	≥ 64 (4/76)
Clarithromycin ³	≤ 8	16	≥ 32
Moxifloxacin	≤ 1	2	≥ 4
Linezolid	≤ 8	16	≥ 32
Minocycline, Doxycycline	< 1	2 – 4	> 8

¹ = CLSI: M24-A2, 2011.

²= Amikacin $>64 \mu\text{g/mL}$ is considered resistant for inhaled treatment; $\geq 64 \mu\text{g/mL}$ is considered resistant using IV treatment. Confirmation of amikacin resistance can be done by ordering mutational sequencing of the 16S rRNA gene (Test Code 13). These MICs are currently being proposed to the CLSI for adoption in the new CLSI M24-A3.

³ = only clarithromycin and amikacin are routinely reported for MAC isolates. Clarithromycin serves as a class drug for all newer macrolides (especially azithromycin).

Please note that the MIC testing of ethambutol (which is considered to be technically difficult and unreliable by many investigators) is currently being discussed by the CLSI and this agent may be removed dependent upon the pending CLSI consensus report.

References

1. Brown BA, Wallace RJ Jr: Broth microdilution MIC test for *Nocardia spp* In: Clinical Microbiology Procedures Handbook, Section 5: Antimicrobial Susceptibility Testing, p. 5.12.1. American Society for Microbiology Book Division, Washington, DC. 1992.
2. Wallace RJ Jr, Nash DR, Steele LC, Steingrube VA: Susceptibility testing of slowly growing mycobacteria utilizing a microdilution MIC method with 7H9 broth . *Journal of Clinical Microbiology* 24:976-981, 1986.
3. Wallace RJ Jr, Tsukamura M, Brown BA, Brown J, Steingrube VA, Zhang Y, Nash DR. Cefotaxime-resistant *Nocardia asteroides* strains are isolates of the controversial species *Nocardia farcinica*. *J. Clin. Microbiol.* 28:2726-2732, December, 1990.
4. McBride ME, Rudolph AH, Tschen JA, Cernoch P, Davis J, Brown BA, Wallace RJ Jr: Diagnostic and therapeutic considerations for cutaneous *Mycobacterium haemophilum* infections. *Arch. Dermatol.* 127:276-277, February, 1991.
5. Wallace RJ Jr., Brown BA, Tsukamura M, Brown JM, Onyi G: Clinical and laboratory features of *Nocardia nova*. *J. Clin. Microbiol.* 29:2407-2411, November, 1991.
6. Burns DN, Wallace RJ Jr, Schultz M, Zhang Y, Zubairi S, Pang Y, Gibert C, Brown BA, Noel E, Gordin FM: Nosocomial outbreak of respiratory tract colonization with *Mycobacterium fortuitum*: Demonstration of the usefulness of pulse-field gel electrophoresis in an epidemiologic investigation. *Am. Rev. Respir. Dis.* 144:1153-1159, 1991.
7. Shapiro CL, Haft RF, Gantz NM, Doern GV, Christenson JC, O'Brien R, Overall JC, Brown BA, Wallace RJ Jr: *Tsukamurella paurometabolum*: A novel pathogen causing catheter-related bacteremia in patients with cancer. *Clin. Infect. Dis.* 14:200-203, January, 1992.
8. Brown BA, Swenson, JM, Wallace RJ Jr: Broth microdilution MIC Test for rapidly growing mycobacteria. In: Clinical Microbiology Procedures Handbook, Sect. 5: Antimicrobial Susceptibility Testing, p. 5.11.1, American Society for MicrobiologIn: Clinical Microbiology Procedures Handbook, Sect. 5: Antimicrobial Susceptibility Testing, p. 5.11.1, American Society for Microbiology Book Division, Washington, DC, 1992.
9. Brown BA, Wallace RJ Jr, Onyi G: Activities of clarithromycin against eight slowly growing species of nontuberculous mycobacteria, determined by using a broth microdilution MIC system. *Antimicrob. Agents Chemother.* 36:1987-1990, September, 1992.
10. Wallace RJ Jr, Brown BA, Blacklock Z, Ulrich R, Jost K, Brown JM, McNeil MM, Onyi G, Steingrube VA, Gibson JL: New *Nocardia* taxon among isolates of *Nocardia brasiliensis* associated with invasive disease. *J. Clin. Microbiol.* 33:1528-1533, June, 1995.
11. Ruimy R, Riegel P, Carlotti A, Boiron P, Bernardin G, Monteil H, Wallace RJ Jr, Christen R: *Nocardia pseudobrasiliensis* sp. nov., a new species of *Nocardia* which groups bacterial strains previously identified as *Nocardia*

- brasiliensis* and associated with invasive diseases. *Intern. J. Syst. Bacteriol.* 46:259-264, January, 1996.
12. Wilson, RW, Steingrube VA, BA Brown, Blacklock Z, Jost Jr. KC, McNabb A, Colby WD, Biehle JR, Gibson JL, Wallace RJ Jr: Recognition of a *Nocardia transvalensis* complex by resistance to aminoglycosides, including amikacin, and PCR-restriction fragment length polymorphism analysis. *J. Clin. Microbiol.* 35:2235-2242, September 1997.
 13. Wallace RJ Jr., Cook JL, Glassroth J, Griffith De, Olivier KN, Gordin F: Diagnosis and treatment of disease caused by nontuberculous mycobacteria . *American Thoracic Society Statement. American Journal of Respiratory and Critical Care Medicine*, 156:S1-S25, August 1997.
 14. Steingrube VA, Gibson JL, Brown BA, Zhang Y, Wilson Endonuclease analysis of a 65-Kilodalton heat shock protein gene sequence for taxonomic separation of rapidly growing mycobacteria. *Journal of Clinical Microbiology*, 33:149-153, 1995.
 15. Steingrube VA, Brown BA, Gibson JL Wilson RW, Brown J, Blacklock Z, Jost K, Locke S, Ulrich RF, Wallace RJ Jr: DNA amplification and restriction endonuclease analysis for differentiation of 12 species and taxa of *Nocardia*, including recognition of four new taxa within the *Nocardia asteroides* complex. *Journal of Clinical Microbiology*, 33:3096-3101, 1995.
 16. Steingrube VA, Wilson RW, Brown BA, Jost KC Jr., Blacklock Z, Gibson JL, Wallace RJ Jr: Rapid identification of clinically significant species and taxa of aerobic actinomycetes, including *Actinomyces*, *Gordona*, *Nocardia*, *Rhodococcus*, *Streptomyces*, and *Tsukamurella* isolates, by DNA amplification and restriction endonuclease analysis. *Journal of Clinical Microbiology*, 35:817-822, April 1997.
 17. Telenti A, Marchesi F, Balz M, Bally F, Böttger EC, Bodmer, T: Rapid identification of mycobacteria to the species level by polymerase chain reaction and restriction enzyme analysis . *Journal of Clinical Microbiology*, 31:175-178, 1993.
 18. Swenson JM, Wallace RJ Jr, Silcox VA, Thornsberry C: Antimicrobial susceptibility of 5 subgroups of *Mycobacterium fortuitum* and *Mycobacterium chelonae*. *Antimicrobial Agents Chemotherapy* 28:807-811, 1985.
 19. Woods GL, Bergmann JS, Witebsky FG, Fahle GA, Wanger A, Boulet B, Plaunt M, Brown BA, Wallace RJ Jr: Multisite reproducibility of results obtained by the broth microdilution method for susceptibility testing of *Mycobacterium abscessus*, *Mycobacterium chelonae*, and *Mycobacterium fortuitum* . *Journal of Clinical Microbiology*, 37:1676-1682, 1999.
 20. Zhang Y, Rajagopalan M, Brown BA, Wallace RJ Jr: Randomly amplified polymorphic DNA PCR for comparison of *Mycobacterium abscessus* strains from nosocomial outbreaks. *Journal of Clinical Microbiology*, 35:3132-3139, 1997.
 21. Hector JSR, Pang Y, Mazurek GH, Zhang Y, Brown Ba, Wallace RJ Jr: Large restriction fragment patterns of genomic *Mycobacterium fortuitum* DNA as strain-specific markers and their use in epidemiologic

- investigation of four nosocomial outbreaks. *Journal of Clinical Microbiology*, 30:1250-1255, 1992.
22. Wallace RJ Jr, Zhang Y, Brown BA, Fraser V, Mazurek GH, Maloney S: DNA large restriction fragment patterns of sporadic and epidemic nosocomial strains of *Mycobacterium chelonae* and *Mycobacterium abscessus*. *Journal of Clinical Microbiology*, 31:2697-2701, 1993.
 23. Villanueva A, Calderon RV, Vargas BA Ruiz F, Aguero S, Zhang Y, Brown BA, Wallace RJ Jr: Report on an outbreak of postinjection abscesses due to *Mycobacterium abscessus*, including management with surgery and clarithromycin therapy and comparison of strains by random amplified polymorphic DNA polymerase chain reaction. *Clinical Infectious Diseases*, 24:1147-1153, 1997.
 24. Lai KK, Brown BA, Westerling JA, Fontecchio SA, Zhang Y, Wallace RJ Jr: Long-term laboratory contamination by *Mycobacterium abscessus* resulting in two pseudo-outbreaks: Recognition with use of random amplified polymorphic DNA (RAPD) polymerase chain reaction. *Clinical Infectious Diseases*, 27:169-175, 1998.
 25. Meier A, Heifets L, Wallace RJ Jr, Zhang Y, Brown BA, Sander P, Böttger EC. Molecular mechanisms of clarithromycin resistance in *Mycobacterium avium*: observation of multiple 23S rDNA mutations in a clonal population. *The Journal of Infectious Diseases*, 174:354-60, 1996.
 26. Brown BA, Lopes JO, Wilson RW, Costa JM, de Vargas AC, Alves SH, Klock C, Onyi GO, Wallace RJ Jr: Disseminated *Nocardia pseudobrasiliensis* infection in a patient with AIDS in Brazil. *Clin. Infect. Dis.* 28:144-145, January 1999.
 27. Brown BA, Springer B, Steingrube VA, Wilson RW, Pfyffer GE, Garcia MJ, Menendez MC, Rodriguez-Salgado B, Jost KC Jr., Chiu SH, Onyi GO, Böttger EC, Wallace RJ Jr: *Mycobacterium wolinskyi* sp. nov. and *Mycobacterium goodii* sp. nov., two new rapidly growing species related to *Mycobacterium smegmatis* and associated with human wound infections: a cooperative study from the International Working Group on Mycobacterial Taxonomy. *Intern. J. Syst. Bacteriol.* 49:1493-1511, 1999.
 28. Wallace RJ Jr, Brown-Elliott BA, Ward SC, Crist CJ, Mann LB, Wilson RW: Activities of linezolid against rapidly growing mycobacteria. *Antimicrob. Agents Chemother.* 45:764-767, March 2001.
 29. Brown-Elliott BA, Ward SC, Crist CJ, Mann LB, Wilson RW, Wallace RJ Jr: In Vitro Activities of linezolid against multiple *Nocardia* species. *Antimicrob. Agents Chemother.* 45:1295-1297, April 2001.
 30. Wilson RW, Steingrube VA, Böttger EC, Springer B, Brown-Elliott BA, Vincent V, Jost KC Jr., Zhang Y, Garcia MJ, Chiu SH, Onyi GO, Rossmore H, Nash DR, Wallace RJ Jr: *Mycobacterium immunogenum* sp. nov., a novel species related to *Mycobacterium abscessus* and associated with clinical disease, pseudo-outbreaks and contaminated

- metalworking fluids: an international cooperative study on mycobacterial taxonomy. *Int. J. Syst. Evol. Microbiol.* 51:1751-1764, 2001.
31. Wallace RJ Jr, Brown-Elliott BA, Hall L, Roberts G, Wilson RW, Mann LB, Crist CJ, Chiu SH, Dunlap R, Garcia MJ, Bagwell JT, Jost KC Jr: Clinical and laboratory features of *Mycobacterium mageritense*. *J. Clin. Microbiol.* 40:2930-2935, August 2002.
 32. Wallace RJ Jr, Brown-Elliott BA, Crist CJ, Mann L, Wilson RW: Comparison of the in vitro activity of the glycolcycline tigecycline (formerly GAR-936) with those of tetracycline, minocycline, and doxycycline against isolates of nontuberculous mycobacteria. *Antimicrob. Agents Chemother.* 46: 3164-3167, October 2002.
 33. Brown-Elliott BA, Wallace RJ Jr: Clinical and taxonomic status of pathogenic nonpigmented or late-pigmenting rapidly growing mycobacteria. *Clin. Microbiol. Rev.* 15:716-746, October 2002.
 34. Woods GL, Williams-Bouyer N, Wallace RJ Jr, Brown-Elliott BA, Witebsky FG, Conville PS, Plaunt M, Hall G, Aralar P, Inderlied C: Multisite reproducibility of results obtained by two broth dilution methods for susceptibility testing of *Mycobacterium avium* complex. *J. Clin. Microbiol.* 41:627-631, February 2003.
 35. Schinsky MF, Morey RE, Steigerwalt AG, Douglas MP, Wilson RW, Floyd MM, Butler WR, Daneshvar MI, Brown-Elliott BA, Wallace RJ Jr, McNeil MM, Brenner DJ, Brown JM: Taxonomic variation in the *Mycobacterium fortuitum* third biovariant complex: description of *Mycobacterium boenickei* sp. nov., *Mycobacterium houstonense* sp. nov., *Mycobacterium neworleansense* sp. nov. and *Mycobacterium brisbanense* sp. nov. and recognition of *Mycobacterium porcinum* from human clinical isolates. *Int. J. Syst. Evol. Microbiol.* 54:1653-1667, March 2004.
 36. Brown-Elliott BA, Brown JM, Conville PS, Wallace RJ Jr: Clinical and laboratory features of the *Nocardia* spp. based on current molecular taxonomy. *Clin. Microbiol. Rev.* 19:259-282, 2006.
 37. Brown-Elliott, BA, Cohen, S, Wallace RJ Jr.: Susceptibility testing of mycobacteria. In: Antimicrobial Susceptibility Testing Protocols, R Schwalbe, L Steele-Moore, AC Goodwin, eds. CRC Press, Boca Raton, FL. 2007; 11:243-274.
 38. Forbes BA, Banaiee N, Beavis KG, Brown-Elliott BA, Della Latta P, Elliott LB, Hall GS, Hanna B, Perkins MD, Siddiqi SH, Wallace RJ Jr, Warren NG: Laboratory detection and identification of mycobacteria; Approved guideline. *Clinical and Laboratory Standards Institute*, M48-A, 2008.
 39. Petti CA, Bosshard PP, Brandt ME, Clarridge JE III, Feldblyum TV, Foxall P, Furtado MR, Pace N, and Procop G. Interpretive Criteria for Identification of Bacteria and Fungi by DNA Target Sequencing; Approved Guideline. *Clinical and Laboratory Standards Institute*, MM18-A, 2008.
 40. Nash, K. A., B. A. Brown-Elliott, and R. J. Wallace Jr. A novel gene, *erm(41)*, confers inducible macrolide resistance to clinical isolates of

- Mycobacterium abscessus* but is absent from *Mycobacterium chelonae*. *Antimicrob. Agents Chemother.* 2009; 53:1367-1376.
41. Brown-Elliott BA, Wallace RJ Jr: Rapidly growing mycobacteria. In: Tuberculosis & Nontuberculous Mycobacterial Infections, 6th edition, D Schlossberg, ed. McGraw-Hill Companies, Inc., 2010, Accepted.
 42. Woods GL, Brown-Elliott BA, Conville PS, Desmond EP, Hall GS, Lin G, Pfyffer GE, Ridderhof JC, Siddiqi SH, Wallace RJ Jr, Warren NG, Witebsky FG: Susceptibility testing of mycobacteria, nocardiae, and other aerobic actinomycetes; Approved Standard - Second Edition. CLSI, 2011 document M24-A2.
 43. Richter E, Brown-Elliott BA, Wallace RJ Jr.: *Mycobacterium*: Laboratory characteristics of slowly growing mycobacteria. In: Manual of Clinical Microbiology, 9th Edition, P.R. Murray, ed. ASM Press, Washington, D.C. 2011, Vol. 1, pp. 503-524.
 44. Brown-Elliott BA, Wallace RJ Jr.: *Mycobacterium*: Clinical and Laboratory Characteristics of rapidly growing mycobacteria. In: Manual of Clinical Microbiology, 10th Edition, P.R. Murray, ed. ASM Press, Washington, D.C. 2011, Vol. 1, pp. 525-538.
 45. Brown-Elliott BA, Crist CJ, Mann LB, Wilson RW, Wallace RJ Jr: vitro activity of linezolid against slowly growing nontuberculous mycobacteria. *Antimicrob. Agents Chemother.* 47:1736-1738, May 2003.
 46. Schinsky MF, Morey RE, Steigerwalt AG, Douglas MP, Wilson RW, Floyd MM, Butler WR, Daneshvar MI, Brown-Elliott BA, Wallace RJ Jr, McNeil MM, Brenner DJ, Brown JM: Taxonomic variation in the *Mycobacterium fortuitum* third biovariant complex: description of *Mycobacterium boenickei* sp. nov., *Mycobacterium houstonense* sp. nov., *Mycobacterium new orleansense* sp. nov and *Mycobacterium brisbanense* sp. nov. and recognition of *Mycobacterium porcinum* from human clinical isolates. *Int. J. Syst. Evolution. Microbiol.* 54:1653-1667, 2004.
 47. Brown-Elliott BA, Wallace RJ Jr: Rapidly growing mycobacteria. In: Tuberculosis & Nontuberculous Mycobacterial Infections, 6th edition, D Schlossberg, ed. McGraw-Hill Companies, Inc. 37:565*577, 2011.
 48. Griffith DE, Aksamit T, Brown-Elliott BA, Catanzaro A, Daley C, Gordin, Holland SM, Horsburgh R, Huitt G, Iademarco MF, Iseman M, Olivier K, Ruoss S, von Reyn FC, Wallace RJ Jr., Winthrop K. An Official ATS/IDSA Statement: Diagnosis, Treatment, and Prevention of Nontuberculous Mycobacterial Diseases. *Am. Respir. Crit. Care Med*, 175:367-416, 2007.
 49. Nash KA, Brown-Elliott, BA, Wallace JR Jr. A Novel Gene *erm*(41), Confers inducible macrolide resistance to clinical isolates of *Mycobacterium abscessus* but is absent from *Mycobacterium chelonae*. *Antimicrob. Agents Chemother*, 53:1367-1376, 2009.
 50. Brown-Elliott BA, Nash KA, Wallace RJ Jr. Antimicrobial susceptibility testing, drug resistance mechanisms, and therapy for infections with nontuberculous mycobacteria. *Clin. Microbiol. Rev.* 25:545-582, 2012.

51. Brown-Elliott BA, Iakhiaeva E, Griffith DE, Woods G, Stout JE, Wolfe CR, Turenne C, Wallace RJ Jr.: *In vitro* activity of amikacin against isolates of *Mycobacterium avium* complex with proposed MIC breakpoints and finding of a 16S rRNA gene mutation in treated isolates. *J. Clin. Microbiol.* 51:3389-3394, 2013.
52. Macheras E, Roux A-L, Ripoll F, Sivadon-Tardy V, Gutierrez C, Gaillard J-L, Heym B: Inaccuracy of single-target sequencing for discriminating species of the *Mycobacterium abscessus* group. *J. Clin. Microbiol.* 47:2596-2600, 2009.
53. Macheras E, Roux A-L, Bastian S, Leão SC, Palaci M, Sivadon-Tardy V, Gutierrez C, Richter E, Rüscher-Gerdes S, Pfyffer G, Bodmer T, Cambau E, Gaillard J-L, Heym B: Multilocus sequence analysis and *rpoB* sequencing of *Mycobacterium abscessus* (sensu lato) strains. *J. Clin. Microbiol.* 49:491-499, 2011.
54. Tettelin H, Davidson RM, Agrawal S, Aitken ML, Shallom S, Hasan NA, Strong M, de Moura VCN, De Groot MA, Duarte RS, Hine E, Parankush S, Su Q, Daugherty SC, Fraser CM, Brown-Elliott BA, Wallace RJ, Holland SM, Sampaio EP, Olivier KN, Jackson M, Zelazny AM: High-level relatedness among *Mycobacterium abscessus* subsp. *massiliense* strains from widely separated outbreaks. Current Opinion, Centers for Disease Control and Prevention, ISSN: 1080-6059, eid/article/20/3/13-1106 article.htm#comment. 20(3), March 2014.
55. Adékambi T, Colson P, Drancourt M: *rpoB*-based identification of nonpigmented and late-pigmenting rapidly growing mycobacteria. *J. Clin. Microbiol.* 41:5699-5708, 2003.
56. Iakhiaeva E, McNulty S, Brown-Elliott BA, Falkinham III JO, Williams MD, Vasireddy R, Wilson RW, Turenne C, Wallace Jr. RJ: Mycobacterial interspersed repetitive-unit-variable-number tandem-repeat (MIRU-VNTR) genotyping of *Mycobacterium intracellulare* for strain comparison with establishment of a PCR-based database. *J. Clin. Microbiol.* 51:409-416, 2013.
57. Conville PS, Zelazny AM, Witebsky FG: Analysis of *secA1* gene sequences for identification of *Nocardia* species. *J. Clin. Microbiol.* 2006; 44:2760-2766.
58. Brown-Elliott BA, Hanson K, Vasireddy S, Iakhiaeva E, Nash KA, Vasireddy R, Parodi N, Smith T, Gee M, Strong A, Baker A, Cohen S, Muir H, Slechta ES, Wallace Jr. RJ: Absence of a function *erm* gene in isolates of *Mycobacterium immunogenium* and the *Mycobacterium mucogenicum* group, based on *in vitro* clarithromycin susceptibility. *J. Clin. Microbiol.* 53:875-878, 2015.
59. Brown-Elliott BA, Vasireddy S, Vasireddy R, Iakhiaeva E, Howard ST, Nash K, Parodi N, Strong A, Gee M, Smith T, Wallace Jr. RJ: Utility of sequencing the *erm(41)* gene in isolates of *Mycobacterium abscessus*

- subsp. *abscessus* with low and intermediate clarithromycin MICs. *J. Clin. Microbiol.* 53:1211-1215, 2015.
60. Nash KA, Andini N, Zhang Y, Brown-Elliott BA, Wallace Jr., RJ: Intrinsic macrolide resistance in rapidly growing mycobacteria. *Antimicrob. Agents Chemother.* 50:3476-3478, 2006.
61. Nash KA, Zhang Y, Brown-Elliott BA, Wallace RJ: Molecular basis of intrinsic macrolide resistance in clinical isolates of *Mycobacterium fortuitum*. *J. Antimicrob. Chemother.* 55:170-177, 2005.
62. Iakhiava E, Howard ST, Brown-Elliott BA, McNulty S, Newman KL, Falkinham III JO, Williams M, Kwait R, Lande L, Vasireddy R, Turenne C, Wallace Jr., RJ: Variable-number tandem-repeat analysis of respiratory and household water biofilm isolates of "*Mycobacterium avium* subsp. *hominissuis*" with establishment of a PCR database. *J. Clin. Microbiol.* 54:891-901.